

Materials List for:

The Modified Hole Board – Measuring Behavior, Cognition and Social Interaction in Mice and Rats

Maaïke Labots^{1,2}, Hein A. Van Lith^{1,2}, Frauke Ohl^{1,2}, Saskia S. Arndt^{1,2}

¹Department of Animals in Science and Society, Utrecht University

²Brain Center Rudolf Magnus

Correspondence to: Maaïke Labots at M.Labots@uu.nl

URL: <http://www.jove.com/video/52529>

DOI: [doi:10.3791/52529](https://doi.org/10.3791/52529)

Materials

Name	Company	Catalog Number	Comments
Cognitive mHB apparatus	N/A	N/A	The box for cognitive testing is 50x50 cm by inserting a partition made of grey PVC in the standard mHB. A smaller board (35x22x1 cm) with 10 cylinders is placed in the middle of the box. All cylinders are scented with a flavor animals are attracted to (e.g. vanilla) and all are baited with a reward (e.g. a piece of almond) beneath a grid so the animals cannot remove it. Cylinders (often three) are cued with a colored ring (contrasting with the grey PVC) and are baited with a removable reward (e.g. 0.05 g piece of almond).
Vanilla flavor	N/A	N/A	Vanilla flavor dissolved in water (0.02%)
Camera	N/A	N/A	Camera that generates sufficient quality output to rescore the behavior from video.
Stage light	N/A	N/A	A theater light that allows to adjust the light intensity and the dimensions of the surface that is lit.
Behavioral scoring software	N/A	N/A	Behavioral scoring software.